Inhaler Assessment Kit

CE 0120

INTERNATIONAL

Clement Clarke

A HAAG-STREIT COMPANY

© Copyright 1998
Clement Clarke International Ltd.

In-check Inhaler Assessment Kit Instruction Booklet – Issue 3 3/99	Ref. 1902042
In-check Inhaler Assessment Kit	Ref. 3109760
In-check Peak Inspiratory Flow Meter	Ref. 3109750
Disposable Cardboard Mouthpieces – Adult, 5 x 20	Ref. 3125030
Disposable Cardboard Mouthpieces – Adult, 500	Ref. 3122003
Universal Plastic Mouthpieces – Sterilisable, 15's	Ref. 3103095

Clement Clarke International Limited
Edinburgh Way
Harlow
Essex
CM20 2TT
U.K.

Tel. +44 (0) 1279 414969
Fax. + 44 (0) 1279 635232
email: info@resp.clement-clarke.com
www: www.inspiratory.com

Contents

English	G	Instructions for Use	Page 2
Français	F	Mode d'emploi	Page 4
Deutsch	D	Gebrauchsanleitung	Page 6
Italiano	I	Istruzioni per l'uso	Page 8
Nederlands	O	Gebruiksaanwijzing	Page 10
Español	E	Instrucciones de Uso	Page 12
Português	P	Instruções de Utilização	Page 14
Íslenska	g	Notkunarreglur	Page 16
Ελληνικά	K	Οδηγίες Χρήσης	Page 18
Suomi	O	Käyttöohjeet	Page 20
Dansk	C	Brugsvejledning	Page 22
Svenska	S	Instruktionsanvisning	Page 24
Norsk	N	Brukerveiledning	Page 26

CE 0120

Attention – Consult accompanying documents.

Caution: Federal law restricts this device to sale by or on the order of a physician.

Before operating the instrument, please read the instructions supplied carefully, and keep them safe for future reference.

This product complies with the essential requirements of the medical devices directive (93/42/EEC). Compliance has been verified by UK notified body per 0120 (SGS Yarsley).

The In-check Inhaler Assessment Kit has been designed to provide an objective measure of the inspiratory flow achieved through various inhaler devices. The efficacy of dry powder inhalers is dependent on the inspiratory flow rate at which they are used, and breath-actuated metered dose inhalers will only trigger release of medication above a specific flow rate. By checking the patient's flow rate for a range of devices, the health professional can identify which inhalers are suitable for use.

The kit contains a hand-calibrated inspiratory flow meter (the **In-check**), a sterilisable universal mouthpiece, and a series of adapters that reproduce the resistance of popular inhaler devices. The plastic universal mouthpiece may be substituted by cardboard disposable mouthpieces (available separately).

A table of minimum (**MINIMUM FLOW**) and optimum flows (**OPTIMUM FLOW**) for specific inhaler devices will be found on the back cover of this booklet.

IMPORTANT

With any inhalation device, it is important to check for loose foreign objects before the device is used. The transparent material used in the construction of the **In-check** enables the user to make a visual check before inhalation. Patients should be prevented from exhaling through the device prior to use.

To reset the In-check

Hold the instrument vertically with the mouthpiece uppermost, so the rounded end of the meter can be tapped from underneath with the other hand.

A hard tap will dislodge the magnetic resetting weight, which will return the cursor to a start position. When this has happened, invert the meter through 180 degrees to return the magnetic weight to its resting position.

How to use the In-check to assess inspiratory flow

1. Reset the **In-check**.
2. Attach an adapter and clean mouthpiece.
3. Ask the patient to exhale fully.
4. Seal lips around the mouthpiece.
5. According to the adapter chosen, instruct the patient to inhale appropriately – as recommended by the relevant inhaler manufacturer (some manufacturers recommend “slow and controlled”, others “fast and strong”).
6. Record the inspiratory flow from the position of the red cursor against the scale. Reset, and repeat two more times.
7. Compare values achieved with target flows for that device – to operate an inhaler device effectively, the patient should be able to achieve the minimum (**MINIMUM FLOW**) value.

For some devices, the best drug deposition is achieved when the patient can achieve the optimal (**OPTIMAL FLOW**) value.
8. If the minimum (**MINIMUM FLOW**) value cannot be achieved, then an alternative type of inhaler device should be assessed.

Performance Accuracy

Accuracy $\pm 10\%$ or 10 l/min (whichever is greater) and repeatability of ± 5 l/min.

Cleaning your In-check

Immerse **In-check** in warm (but not hot) mild detergent solution for 2-3 minutes (maximum 5 minutes). Agitate the meter to ensure thorough cleaning.

Rinse in warm water and shake to remove any excess water. It is important to rinse thoroughly to prevent salt spots appearing on the inside of the body and the spindle.

Allow to dry thoroughly before using again.

The **In-check** restrictive adapter should also be cleaned using the procedure above, ensuring that the central orifice is completely clean.

The expected life of the **In-check** meter, in normal use, is two years.

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References:

1. Pedersen S, Hansen O R, Fuglsang G. Influence of inspiratory flow rate upon the effect of a Turbuhaler™. *Arch Dis Child* 1990; **65**: 308-319
2. Dolovich M B, Vanzeileghem M, Hidinger K-G, Newhouse M. Influence of inspiratory flow rate (VI) on the response to terbutaline (T) via the Turbuhaler™. *Allergy*, 1989; **44**: 220-225
3. Persson G, Gruvstad E, Stahl E. A new multiple dose powder inhaler, (Turbuhaler™), compared with a pressurised inhaler in a study of terbutaline in asthmatics. *Eur Respir J*, 1988; **1**: 681-684
4. Borgstrom L, Bondesson E, Moren F et al. Lung deposition of budesonide inhaled via Turbuhaler™: a comparison with terbutaline sulphate in normal subjects. *Eur Respir J* 1994; **7**: 69-73
5. Johnsen C R, Weeke E R. Turbuhaler™: a new device for dry powder terbutaline inhalation. *Allergy*, 1988; **43**: 392-395.
6. Borgstrom L, Derom E, Stahl E, Wahlin-Bull E, Pauwels R. The inhalation device influences lung deposition and bronchodilating effect of terbutaline. *Am J Respir Crit Care Med* 1996; **153**: 1436-1440
7. Hill L S, Slater A L. A comparison of the performance of two modern multidose dry powder asthma inhalers. *Resp Med* 1998; **92**: 105-110
8. Nielsen K G, Auk I L, Bojsen K et al. Clinical effect of Diskus™ dry-powder inhaler at low and high inspiratory flow rates in asthmatic children. *Eur Resp J* 1998; **11**: 350-4
9. Malton A, Sumby B S, Dandiker Y. A comparison of in-vitro drug delivery from salbutamol Diskus™ and terbutaline Turbuhaler™ inhalers. *J. Pharm. Med* 1996; **6**: 35-48
10. Malton A, Sumby B S, Smith I J. A comparison of in-vitro drug delivery from two multidose powder inhalation devices. *Eur J Clin Res* 1995; **7**: 177-193
11. Prime D, Sharma R K, Slater A L, Petchley L, Parkes P A. In-vitro comparison of the dose delivery characteristics of two multi-dose powder inhalers delivering long-acting B₂-agonists. *Eur Resp J* 1997; **10** (suppl 25): 258s
12. Personal communication. Re: Data on file. Novartis Pharmaceuticals UK Ltd. August 1998
13. Nielsen KG, Skov M, Klug B, et al. Flow dependent effect of formoterol dry-powder inhaled from the Aerolizer®. *Eur Respir J* 1997; **10**: 2105-2109
14. Baum E A, Bryant A M. The development and laboratory testing of a novel breath-actuated pressurized inhaler. *J Aerosol Med* 1988; **1**: 219-220
15. Newman S P, Weisz A W B, Talaei N, Clarke S W. Improvement of drug delivery with a breath actuated pressurised aerosol for patients with poor inhaler technique. *Thorax* 1991; **46** (10): 712-716
16. Ross D L, Schultz R K. Effect of inhalation flow rate on the dosing characteristics of dry powder inhaler (DPI) and metered dose inhaler (MDI) products. *Journal of Aerosol Medicine* 1996; **9**, (2): 215-226
17. Newman S P, Pavia D, Garland N, Clarke S W. Effects of various inhalation modes on the deposition of radioactive pressurized aerosols. *Europ J Respir Disease* 1982; **63** (SUPPL 119): 57-65
18. Hardy J G, Jasuja A K, Frier M, Perkins A C. A small volume spacer for use with a breath operated pressurised metered dose inhaler. *Int. J. Pharmaceutics* 142 (1996) 129-133

Table of Minimum and Optimum Inspiratory Flows

DPI Devices	Adapter ID	Minimum Flow Manufacturer's recommended minimum flow rate (L/min)	Optimum Flow Flow rate at which optimum dose observed (L/min)	Variation in dose over range
Turbulent flow inhaler	T	30 ref 1,2	60 ref 2, 3, 4, 5	High ref 6,7
Multiple-dose powder inhaler	A/A/D	30 ref 7,8,9,10,11	30 ref 7,8,9,10,11	Low ref 7,8,9,10,11
Aerosol capsule inhaler	F	60 ref 12	120 ref 13	High ref 13
Breath actuated MDIs				
Auto inhaler	A/A/D	30 ref 14,15	30 ref 16,17	v. Low ref 16
Breath-Actuated inhaler	E/S	20 ref 18	20 ref 18	N/A

Date of preparation: March 1999.

Please contact the device manufacturer for additional information.