
easy care

CE

de	5
en	9
bg	13
cs	17
da	21
el	25
es	29
et	33
fi	37
fr	41
hu	45
it	49
lt	53
lv	57
nl	61
no	65
pl	69
pt	73
ro	77
ru	81
sk	85
sl	89
sv	93
tr	97
he	101

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
easy care Cleaning Aid für eFlow® Aerosol Head

Part no. 078G6100

Read these instructions **carefully and completely** before using for the first time. Keep these instructions in a safe place.


1 Product description

Hinged cleaning aid for cleaning the eFlow® aerosol head with a built-in funnel opening and receptacle for filling 5 ml isotonic sodium chloride solution (0.9%). The nebuliser connection cord and the control unit of the eFlow® appliance are also required for operation. Information on connecting the control unit to the power supply and on operation can be found in the Instructions for Use of the eFlow® appliance.

2 Intended use

The easycare cleaning aid can contribute to a shortening of the nebulisation time if the performance of the eFlow® aerosol head is deteriorating.

As a result of easycare cleaning, the pores in the membrane of the eFlow® aerosol head are mechanically rinsed free by means of backwashing in the opposite direction.

3 Important information for use

Cleaning with the easycare cleaning aid is recommended if the nebulisation time for the same fill volume of a given substance significantly increases. Additional information can be found in the current version of the corresponding eFlow® Instructions for Use in the Chapter "Troubleshooting", key word "prolonged nebulisation times" ([www.pari.de/Products/Pari Inhalation devices/eFlow](http://www.pari.de/Products/Pari%20Inhalation%20devices/eFlow)).

Example eFlow® rapid

- Measure the nebulisation time with 2.5 ml isotonic sodium chloride solution (0.9%) without inhaling.
- If the nebulisation time takes longer than **4 minutes**, cleaning with the easycare cleaning aid is recommended.

Frequency of use

If needed, no more than 2 times weekly after inhalation and after cleaning the aerosol head (see cleaning steps below).

Note:
More frequent use is **not** beneficial, but can lead to premature material wear.

During its entire period of use, an aerosol head should not be backwashed more than 25 times with the easycare cleaning aid.

The easycare cleaning aid does not extend the life cycle of the aerosol head.

Backwash liquid

5 ml isotonic sodium chloride solution (0.9%)

Cleaning steps

The easycare backwashing procedure is performed as an **additional** measure for the normal cleaning of the aerosol head.

Please ensure that the following cleaning steps are observed during easycare backwashing:

1. Cleaning nebuliser parts and aerosol head^{*)}
2. **Backwashing the aerosol head with easycare**
3. Disinfecting/sterilising nebuliser parts and aerosol head^{*)}

^{*)} Points 1 and 3 see Chapter "Hygienic Re-use" in the eFlow® Instructions for Use


4 Preparation and assembly

4.1 Open the easycare cleaning aid


- To open press lightly on the lateral gripping surfaces.
- Separate the parts.

4.2 Install the aerosol head and close the easycare reservoir


- The eFlow[®] aerosol head must be felt to snap into place upon insertion. Please take care to ensure proper orientation.
- When closing ensure that the two lateral engaging hooks are felt to snap into place.

4.3 Connect the easycare cleaning aid to the eFlow[®] control unit

- Attach the nebuliser connector of the eFlow[®] nebuliser cable to the hinge of the easycare reservoir and swivel it downwards, so that the metal contacts engage.
- Connect the other end of the cable to the eFlow[®] control unit.

4.4 Filling

- Fill approx. 5 ml isotonic sodium chloride solution (0.9%) into the easycare funnel opening.

5 Backwash procedure

5.1 Operation

- Start the backwash by pressing the start button on the eFlow® control unit.
The processing time for 5 ml is at least 10 min.

5.2 End

After rinsing, the eFlow® control unit automatically switches off. If liquid is still present in the funnel, please restart.

Note:

i At the end of the rinsing process, it is possible for mist to escape from the easycare filling funnel. This is normal. The eFlow® control unit then can also be switched off manually by pressing the ON/OFF button.

6 Disassemble the easycare cleaning aid

- Detach the eFlow® connection cord from the easycare reservoir and from the eFlow® control unit.
- Open the easycare reservoir and remove the eFlow® aerosol head by pressing on the two lateral engaging hooks (disinfection or sterilisation of the aerosol head see Chapter “Hygienic Re-use” in the eFlow® Instructions for Use).
- Discard the rinsing solution from the easycare reservoir.

7 Clean the easycare cleaning aid

After each use, clean with warm tap water and dishwashing liquid. Then dry thoroughly and store opened and turned over on a clean surface.

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PARI Pharma GmbH