

Normalwerte für FEV1 nach EGKS, übernommen von der ERS 1994

		Größe	135	140	145	150	155	160	165	170	175	180	185	190	195
	Alter														
Buben	5	1.93	2.14	2.36	2.6	2.85	3.11	3.39	3.69	4,01	4.33	4.68	5.04	5.42	
	10	1.94	2.14	2.36	2.6	2.85	3.11	3.39	3.69	4,01	4.33	4.68	5.04	5.42	
	15	1.95	2.14	2.36	2.6	2.85	3.11	3.39	3.69	4,01	4.33	4.68	5.04	5.42	
Männer	20	2.58	2.8	3.01	3.23	3.44	3.66	3.87	4.09	4.3	4.52	4.73	4.95	5.16	
	25	2.58	2.8	3.01	3.23	3.44	3.66	3.87	4.09	4.3	4.52	4.73	4.95	5.17	
	30	2.44	2.65	2.87	3.08	3.3	3.51	3.73	3.94	4.16	4.37	4.59	4.8	5.02	
	35	2.29	2.51	2.72	2.94	3.15	3.37	3.58	3.8	4.01	4.23	4.44	4.66	4.87	
	40	2.15	2.36	2.58	2.79	3.01	3.22	3.44	3.65	3.87	4.08	4.3	4.51	4.73	
	45	2	2.22	2.43	2.65	2.86	3.08	3.29	3.51	3.72	3.94	4.15	4.37	4.58	
	50	1.86	2.07	2.29	2.5	2.72	2.93	3.15	3.36	3.58	3.79	4.01	4.22	4.44	
	55	1.71	1.93	2.14	2.36	2.57	2.79	3	3.22	3.43	3.65	3.86	4.08	4.29	
	60	1.57	1.78	2	2.21	2.43	2.64	2.86	3.07	3.29	3.5	3.72	3.93	4.15	
	65	1.42	1.64	1.85	2.07	2.28	2.5	2.71	2.93	3.14	3.36	3.57	3.79	4	
	70	1.28	1.49	1.71	1.92	2.14	2.35	2.57	2.78	3	3.21	3.43	3.64	3.86	
	75	1.13	1.35	1.56	1.78	1.99	2.21	2.42	2.64	2.85	3.07	3.28	3.5	3.71	
	80	0.99	1.2	1.42	1.63	1.85	2.06	2.28	2.49	2.71	2.92	3.14	3.35	3.57	
Mädchen	5	1.93	2.14	2.36	2.6	2.85	3.11	3.39	3.69	4,01	4.33	4.68	5.04	5.42	
	10	1.94	2.14	2.36	2.6	2.85	3.11	3.39	3.69	4,01	4.33	4.68	5.04	5.42	
	15	1.95	2.14	2.36	2.6	2.85	3.11	3.39	3.69	4,01	4.33	4.68	5.04	5.42	
Frauen	20	2.1	2.3	2.5	2.7	2.89	3.09	3.29	3.49	3.68	3.88	4.08	4.28	4.47	
	25	2.1	2.3	2.5	2.7	2.89	3.09	3.29	3.49	3.68	3.88	4.08	4.28	4.47	
	30	1.98	2.18	2.37	2.57	2.77	2.97	3.16	3.36	3.56	3.76	3.95	4.15	4.35	
	35	1.85	2.05	2.25	2.45	2.64	2.84	3.04	3.24	3.43	3.63	3.83	4.03	4.22	
	40	1.73	1.93	2.12	2.32	2.52	2.72	2.91	3.11	3.31	3.51	3.7	3.9	4.1	
	45	1.6	1.8	2	2.2	2.39	2.59	2.79	2.99	3.18	3.38	3.58	3.78	3.97	
	50	1.48	1.68	1.87	2.07	2.27	2.47	2.66	2.86	3.06	3.26	3.45	3.65	3.85	
	55	1.35	1.55	1.75	1.95	2.14	2.34	2.54	2.74	2.93	3.13	3.33	3.53	3.72	
	60	1.23	1.43	1.62	1.82	2.02	2.22	2.41	2.61	2.81	3.01	3.2	3.4	3.6	
	65	1.1	1.3	1.5	1.7	1.89	2.09	2.29	2.49	2.68	2.88	3.08	3.28	3.47	
	70	0.98	1.18	1.37	1.57	1.77	1.97	2.16	2.36	2.56	2.76	2.95	3.15	3.35	
	75	0.85	1.05	1.25	1.45	1.64	1.84	2.04	2.24	2.43	2.63	2.834	3.03	3.22	
	80	0.73	0.93	1.12	1.32	1.52	1.72	1.91	2.11	2.31	2.51	2.7	2.9	3.1	

Ferdinand Menzl
Medizintechnik GmbH

Ferdinand Menzl Medizintechnik GmbH

A-1220 Wien, Donaufelderstraße 199

FON: +43 (1) 2558960-0

FAX: +43 (1) 2558960-0

WEB: www.menzl.com